# EMSCULPT neo<sup>®</sup>

The **ONE** in Body Shaping

MORE FOR LESS

EMSCULPT<sup>neo</sup><sup>®</sup>

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\*Abdominal toning and reduction of subcutaneous fat with combination of novel radiofrequency treatment and HIFEM procedure - MRI scan study. Jacob C et al. Presented at ASDS 2020 Virtual Meeting. Radiofrequency heating and HIFEM delivered simultaneously - the first sham-controlled randomized trial. Katz B et al. Presented at ASDS 2020 Virtual Meeting. Ultrasound evaluation of the simultaneous RF and HIFEM treatments on human fat tissue. Denkova R. Source: U.S. Food and Drug Administration. 510(k) Premarket Notification: K192224. Published online on December 5, 2019. TÊ

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
<sup>1</sup>Hyperthermic Injury to Adipocyte Cells by Selective Heating of Subcutaneous Fat With a Novel Radiofrequency Device: Feasibility Studies. Franco W. et al. Lasers Surg Med 42:361-370 (2010). <sup>2</sup>Weiss RA, Bernardy J, Tichy F. Radiofrequency Treatment Used in Combination with HIFEM therapy: Histological Analysis including Scanning Electron Microscopy of Adipocytes. In: Phoenix, Arizona, USA: ASLMS; 2020.

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

Halaas Y, Bernardy J, Ondrackova P, Dinev I. The skeletal muscle satellite cell activation by a combination of HIFEM procedure and radiofrequency treatment for body contouring: A first look at the NCAM/CD56 facilitated detection by fluorescent microscopy. In: Phoenix, Arizona, USA: ASLMS; 2020. Mayo Clinic 2019. <https://www.mayoclinic.org/healthy-lifestyle/fitness/in-depth/exercise/art-20045517> ; \*Data on file

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

BEFORE

AFTER **EMSCULPT**<sup>neo</sup>

3 MONTHS AFTER THE LAST TREATMENT, COURTESY OF: BRUCE E. KATZ, M.D.

BEFORE

AFTER **EMSCULPT**<sup>neo</sup>

1 MONTH AFTER THE LAST TREATMENT, COURTESY OF: DAVID KENT, M.D.

MORE FOR LESS

**EMSCULPT**<sup>neo</sup>

BEFORE

AFTER **EMSCULPT**<sup>neo</sup>

3 MONTHS AFTER THE LAST TREATMENT, COURTESY OF: BRUCE E. KATZ, M.D.

BEFORE

AFTER **EMSCULPT**<sup>neo</sup>

1 MONTH AFTER THE LAST TREATMENT, COURTESY OF: DAVID KENT, M.D.

MORE FOR LESS

**EMSCULPT**<sup>neo</sup>

BEFORE

AFTER **EMSCULPT**<sup>neo</sup>

3 MONTHS AFTER THE LAST TREATMENT, COURTESY OF: CAROLYN JACOB, M.D.

BEFORE

AFTER **EMSCULPT**<sup>neo</sup>

1 MONTH AFTER THE LAST TREATMENT, COURTESY OF: DAVID KENT, M.D.

MORE FOR LESS

**EMSCULPT**<sup>neo</sup>

BEFORE

AFTER **EMSCULPT**<sup>neo</sup>


1 MONTH AFTER THE LAST TREATMENT, COURTESY OF: DAVID KENT, M.D.

BEFORE

AFTER **EMSCULPT**<sup>neo</sup>

3 MONTHS AFTER THE LAST TREATMENT, COURTESY OF: DAVID KENT, M.D.

MORE FOR LESS

**EMSCULPT**<sup>neo</sup>

BEFORE

AFTER **EMSCULPT**<sup>neo</sup>

1 MONTH AFTER THE LAST TREATMENT, COURTESY OF: DAVID KENT, M.D.

BEFORE

AFTER **EMSCULPT**<sup>neo</sup>

1 MONTH AFTER THE LAST TREATMENT, COURTESY OF: BRUCE E. KATZ, M.D.

MORE FOR LESS

**EMSCULPT**<sup>neo</sup>

BEFORE

AFTER **EMSCULPT**<sup>neo</sup>

IMMEDIATELY AFTER THE LAST TREATMENT, COURTESY OF: KLAUS HOFFMANN, M.D.

MORE FOR LESS

**EMSCULPT**<sup>neo</sup>

BEFORE

AFTER **EMSCULPT<sup>neo</sup>**

1 MONTH AFTER THE LAST TREATMENT, COURTESY OF: BTL AESTHETICS

MORE FOR LESS

**EMSCULPT<sup>neo</sup>**

BEFORE

AFTER **EMSCULPT**<sup>neo</sup>

IMMEDIATELY AFTER THE LAST TREATMENT, COURTESY OF: MELANIE PALM, M.D.

MORE FOR LESS

**EMSCULPT**<sup>neo</sup>

BEFORE

AFTER **EMSCULPT**<sup>neo</sup>

1 MONTH AFTER THE LAST TREATMENT, COURTESY OF: BTL AESTHETICS

MORE FOR LESS

**EMSCULPT**<sup>neo</sup>

BEFORE

AFTER **EMSCULPT**<sup>neo</sup>


1 MONTH AFTER THE LAST TREATMENT, COURTESY OF: BRIAN KINNEY, M.D.


MORE FOR LESS

**EMSCULPT**<sup>neo</sup>

BEFORE

AFTER **EMSCULPT**<sup>neo</sup>


1 MONTH AFTER THE LAST TREATMENT, COURTESY OF: YAEL HALAAS, M.D.


MORE FOR LESS

**EMSCULPT**<sup>neo</sup>

BEFORE

AFTER **EMSCULPT**<sup>neo</sup>


1 MONTH AFTER THE LAST TREATMENT, COURTESY OF: DIANE DUNCAN, M.D.

MORE FOR LESS

**EMSCULPT**<sup>neo</sup>

BEFORE

AFTER **EMSCULPT**<sup>neo</sup>

1 MONTH AFTER THE LAST TREATMENT, COURTESY OF: DIANE DUNCAN, M.D.

MORE FOR LESS

**EMSCULPT**<sup>neo</sup>

BEFORE

AFTER **EMSCULPT**

2 WEEKS AFTER THE LAST TREATMENT, COURTESY OF: SUH ESTHETICS

BEFORE

AFTER **EMSCULPT**

2 MONTHS AFTER THE LAST TREATMENT, COURTESY OF: WINSTON SALEM DERMATOLOGY & SURGERY CENTER

MORE FOR LESS

**EMSCULPT**<sup>neo</sup>

BEFORE


AFTER **EMSCULPT**


4 WEEKS AFTER 4<sup>th</sup> TREATMENT, COURTESY OF: RAMINDER SALUJA, M.D.

BEFORE

AFTER **EMSCULPT**

4 WEEKS AFTER 4<sup>th</sup> TREATMENT, COURTESY OF: RAMINDER SALUJA, M.D.

MORE FOR LESS

**EMSCULPT**<sup>neo</sup>

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MORE FOR LESS

EMSCULPT<sup>neo</sup>®

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**9**

**New  
Clinical  
Studies**

EMSCULPT NEO showed consistency in eliminating fat and building muscle in all latest clinical studies.

**More  
Than 500  
Patients\***

**30**

**Peer-Reviewed  
HIFEM  
Publications**

More than 30 scientific publications since 2018 make HIFEM the most intensively researched technology used in non-invasive body shaping.

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
EMSCULPT<sup>®</sup>  
neo

# HUMAN HISTOLOGY STUDY


Clinical Evaluation of the BTL-899 Device for Non-invasive Lipolysis Human Fat Tissue.  
Denkova R. 2019.

N = 8 subjects

- Biopsies were taken at baseline, 10 days and 20 days post treatment
- **Apoptotic and lipolytic changes observed at 10 days.** At 20 days intensive fat cell disruption peaked
- Deformed nucleus and pyknotic **nucleus indicating cell death**
- Proved the efficacy of EMSCULPT NEO for non-invasive breakdown of fat

Baseline

20 days after

[https://www.accessdata.fda.gov/cdrh\\_docs/pdf19/K192224.pdf](https://www.accessdata.fda.gov/cdrh_docs/pdf19/K192224.pdf)



MORE FOR LESS

EMSCULPT<sup>®</sup>  
neo

# HUMAN ULTRASOUND STUDY

Effect of the BTL-899 Therapy for Non-invasive Lipolysis and Circumference Reduction of Abdomen. Denkova R. 2018.

N = 42 patients (29 women, 13 men)  
3 treatments (30 min each; 1 per week)

[https://www.accessdata.fda.gov/cdrh\\_docs/pdf19/K192224.pdf](https://www.accessdata.fda.gov/cdrh_docs/pdf19/K192224.pdf)

MORE FOR LESS

EMSCULPT<sup>®</sup> neo

# PORCINE MUSCLE HISTOLOGY



Activation of Skeletal Muscle Satellite Cells by a Device Simultaneously Applying High-Intensity Focused Electromagnetic Technology and Novel RF Technology: Fluorescent Microscopy Facilitated Detection of NCAM/CD56 Y. Halaas MD, D. Duncan MD, MVDr. J. Bernardy Ph.D., MVDr. P. Ondrackova Ph.D., I. Dinev DVM, Ph.D., DSc

- 5 white large pigs, 3x30min treatments on abdomen
- The treatment **increases the satellite cells levels by 30.2%** 2 weeks post-treatment
- Histology showed **hypertrophic changes along with signs of newly formed myofibers**
- The effect is **comparable with 12-16 weeks long exercise programs**

Immunofluorescence images of satellite cell levels (red) at baseline (left) and 2 weeks post-treatment (right)

30.2%

Increase in the satellite cell level

[https://academic.oup.com/asi/advance-article/doi/10.1093/asi/sjab002/6089068?searchresult=](https://academic.oup.com/asi/advance-article/doi/10.1093/asi/sjab002/6089068?searchresult=1)

MORE FOR LESS


EMSCULPT<sup>®</sup>  
neo

# PORCINE FAT HISTOLOGY

Radiofrequency Treatment Used in Combination with HIFEM therapy: Histological Analysis including Scanning Electron Microscopy of Adipocytes. R. A. Weiss MD, FAAD, MVDr J. Bernardy Ph.D., Prof MVDr F. Tichy CSc.


- 7 White Large pigs, 3x30min treatments on abdomen
- Post-treatment the adipocytes were **damaged (apoptosis)** and/or **shrunk (lipolysis)**
- **31.1%** decrease of adipocytes' area at 2 weeks post-treatment
- The fat temperature was maintained **just below 45°C for more than 20 minutes**
- Histology and SEM showed **efficacy** of the procedure for **disruption of adipocytes**

BASELINE


Healthy fat cells with well-defined shape

4 DAYS AFTER

Shrunk adipocytes with noticeable membrane ruptures

2 WEEKS AFTER

Disrupted adipocytes with extrusion of lipid droplets

MORE FOR LESS

EMSCULPT<sup>®</sup>  
neo


# HUMAN FAT HISTOLOGY

**Deletion of adipocytes induced by a novel device simultaneously delivering synchronized radiofrequency and HIFEM: Human histological study.**

David J. Goldberg MD, JD


- **Elimination of adipocytes** and significant reduction in size of fat cells resulted in overall reduction of fat tissue
- Adipocytes were **shrunk** with noticeable **ruptures of cytoplasmic membrane**
- **Control samples** demonstrated **no treatment-related changes** in fat tissue
- **Fat temperature** was maintained between **43-45°C** for the entire treatment

**BASELINE**

Healthy fat cells with well-defined shape

**1 MONTH AFTER**

Disrupted (C) and shrunk fat cells

<https://onlinelibrary.wiley.com/doi/10.1111/jocd.13970>

MORE FOR LESS

EMSCULPT<sup>®</sup>  
neo

# THE MRI STUDY

Abdominal toning and reduction of subcutaneous fat with combination of HIFEM procedure and radiofrequency treatment. Carolyn Jacob MD, David Kent MD. N = 41 patients

## RESULTS AT 3 MONTHS POST-TREATMENT

30.8%

Average fat reduction

26.1%


Average increase in muscle

18.8%


Average reduction in abdominal separation

5.9cm

Average circumferential reduction

Before

3 months after

MORE FOR LESS

EMSCULPT<sup>neo</sup>

# THE ULTRASOUND STUDY

Novel radiofrequency device used in combination with HIFEM procedure for abdominal body shaping: Sham-controlled randomized trial. B.E. Katz MD, R. Weiss MD, J. Samuels MD. N = 72 patients.

## RESULTS AT 3 MONTHS

28.3%


Average fat reduction

24.2%


Average increase in muscle

93.9%

Patient satisfaction

Before

1 month after


MORE FOR LESS

EMSCULPT<sup>®</sup> neo

# DOCUMENTED RESULTS OVER TIME

Data based on results of 3 clinical studies and 125+ patients:

Best results were documented 3 months post-treatment.

Abdominal toning and reduction of subcutaneous fat with combination of novel radiofrequency treatment and HIFEM procedure - MRI scan study. Jacob C. et al. Presented at ASDS 2020 Virtual Meeting.  
Radiofrequency heating and HIFEM delivered simultaneously - the first shamcontrolled randomized trial. Katz B. et al. Presented at ASDS 2020 Virtual Meeting.  
Ultrasound evaluation of the simultaneous RF and HIFEM treatments on human fat tissue. Denkova R. Source: U.S. Food and Drug Administration. 510(k) Premarket Notification: K192224. Published online on December 5, 2019.

MORE FOR LESS

EMSCULPT<sup>®</sup> neo

# ESTUDO LATERAL DA COXA


Simultaneous emission of synchronized radiofrequency and HIFEM energy for treatment of lateral thighs: Interim results of the MRI multicentre study. M. Palm MD, MBA, B.M. Kinney MD, FACS, MSME, Y. Halaas MD, FACS, Richard Goldfarb MD, FACS

- **30 subjects** (29-65 y/o, 19.0-34.5 kg/m<sup>2</sup>)
- Four bilateral treatments combining HIFEM+RF on lateral thighs (“Saddlebags”)
- Interim 1-month MRI data showed **reduction of fat thickness by 1.40 cm**
- Hip **circumference decreased** on average by 3.0 cm


93%

Patient satisfaction

Baseline

1 month

MORE FOR LESS

EMSCULPT<sup>neo</sup>

# THE INNER THIGH STUDY

Efficacy and safety of simultaneous application of HIFEM and synchronized radiofrequency for non-invasive lipolysis in inner thighs: Preliminary data. D. Duncan MD, FACS

- **16 subjects** (24-69 y/o, 21.3-35.0 kg/m<sup>2</sup>)
- Four bilateral treatments with small applicators placed over the inner thigh region
- Preliminary data showed **reduced fat thickness by 0.84 cm** at 1 month and **1.02 cm** at 3 months post treatment
- Thigh **circumference decreased by 1.0 cm** and **1.2 cm** at 1-month and 3-months respectively

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








MORE FOR LESS

EMSCULPT<sup>®</sup> neo


FIGURE 3: (a) T1-weighted transversal MRI image shows a minimal hyperintense area where the rupture in the adductor longus was seen on previous MRI images after the initial trauma. (b) T1-weighted coronal MRI image, additional to (a).

MORE FOR LESS

EMSCULPT<sup>neo</sup>

# BICEPS

PRESET	HIFEM Intensity	Initial RF intensity	Treatment time	Number of treatments	Frequency of treatments	Target skin temperature
<b>HIFEM+RF Gentle</b>	Maximum tolerable	100%. Re-adjusted according to patient's feedback	30 min	4	5-10 days apart	40°C
<b>HIFEM Classic</b>	Máxima Tolerada	-	30 min	4	2-3dias de intervalo	-

MORE FOR LESS

EMSCULPT<sup>neo</sup>

# TRICEPS


PRESET	HIFEM Intensity	Initial RF intensity	Treatment time	Number of treatments	Frequency of treatments	Target skin temperature
<b>HIFEM+RF Gentle</b>	Máxima Tolerada	100%. Reajustada de acordo com o feedback do paciente	30 min	4	5-10dias de intervalo	40°C
<b>HIFEM Classic</b>	Máxima Tolerada	-	30 min	4	2-3dias de intervalo	-

MORE FOR LESS

EMSCULPT<sup>neo</sup>

# PANTORILHA


PRESET	HIFEM Intensity	Initial RF intensity	Treatment time	Number of treatments	Frequency of treatments	Target skin temperature
<b>HIFEM+RF Gentle</b>	Máxima Tolerada	100%. Reajustada de acordo com o feedback do paciente	30 min	4	5-10dias de intervalo	40°C
<b>HIFEM Classic</b>	Máxima Tolerada	-	30 min	4	2-3dias de intervalo	-

MORE FOR LESS

EMSCULPT<sup>neo</sup>

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MORE FOR LESS

EMSCULPT<sup>®</sup>  
neo

**OBRIGADA**

MORE FOR LESS

**EMSCULPT** neo®

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©2012